



Recombinant Inner membrane transport permease yadH (yadH)

Product Code	CSB-CF365329EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFN7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MMHLYWVALKSIWAKEIHRFMRIWVQTLVPPVITMTLYFIIFGNLIGSRIGDMHG FSYMQ FIVPGLIMSVITNAYANVASSFFGAKFQRNIEELLVAPVPTHVIIAGYVGGGVA RGLFV GILVTAISLFFVPFQVHSWVFVALTLVLTAVLFSLAGLLNGVFAKTFDDISLVPTF VLTP LTYLGGVFYSLTLLPPFWQGLSHLNPIVYMISGFRYGFGLGINDVPLVTTFTGVLVV FIVAF YLICWSLIQRGRGLRS
Source	in vitro E.coli expression system
Target Names	yadH
Protein Names	Recommended name: Inner membrane transport permease yadH
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.