



Recombinant Bacillus subtilis Uncharacterized transporter YwcJ (ywcJ)

Product Code	CSB-CF331049BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P39608
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	METQALQKVEQYALKKQNI FASSKIRYVLR SILASIFIGFGITAASKTGSYFFMAD SPFAFPAAAVTFGAAILMIAYGGGDLFTGNTFYFTYTALRKKISWRDTLYLWMS SYAGNLIGAILFAILISATGLFEEPSVHSFLIHLAEHKMEPPASELFFRGMLCNWL VCLAFFIPMSLKGEGAKLFTMMLFVFCFFISGFEHSIANMCTFAISLLIEHPDVT LMGAVRNLIPTLGNLTAGIVMMGWMYYTLNPDQ
Source	in vitro E.coli expression system
Target Names	ywcJ
Protein Names	Recommended name: Uncharacterized transporter ywcJ
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.