



Recombinant Escherichia coli Zinc transporter ZupT (zupT)

Product Code	CSB-CF469408ENW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6I411
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain SE11)
Sequence	MSVPLILTILAGAATFIGAFLGVLGQKPSNRLLAFSLGFAAGIMLLISLMEMLPAA LAAE GMSPVLGYGMFIFGLLGYFGLDRMLPHAHPQDLMQKSVQPLPKSIKRTAILLTL GISLHN FPEGIATFVTASSNLELGFGIALAVALHNIPEGLAVAGPVYAATGSKRTAILWAG ISGLA EILGGVLAWLILGSMISPVVMAAIMAAVAGIMVALSVDELMPLAKEIDPNNNPSY GVLCG MSVMGFSLVLLQTAGIG
Source	in vitro E.coli expression system
Target Names	zupT
Protein Names	Recommended name: Zinc transporter ZupT
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.