



Recombinant Escherichia coli O6:K15:H31 Zinc transporter ZupT (zupT)

Product Code	CSB-CF610653EGY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0TD62
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Sequence	MSVPLILTILAGAATFIGAFLGVLGQKPSNRLLAFSLGFAAGIMLLISLMEMLPAA LAAE GMSPVLGYGMFIFGLLGYFGLDRMLPHAHPQDLMQKSVQPLPKSIKRTAILLTL GISLHN FPEGIATFVTASSNLELGFGIALAVALHNIPEGLAVAGPVYAATGSKRTAILWAG ISGLA EILGGVLAWLILGSMISPVVMAAIMAAVAGIMVALSVDELMPPLAKEIDPNNNPSY GVLCG MSVMGFSLVLLQTAGIG
Source	in vitro E.coli expression system
Target Names	zupT
Protein Names	Recommended name: Zinc transporter ZupT
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.