



Recombinant Human Gamma-secretase subunit APH-1B (APH1B)

Product Code	CSB-CF819900HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8WW43
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTAAVFFGCAFIAGFPALALYVFTIATEPLRIIFLIAGAFFWLVSLLISSLVWFMAR VII DNKDGPTQKYLLIFGAFVSVYIQEMFRFAYYKLLKKASEGLKSINPGETAPSMR LLAYVS GLGFGIMSGVFSFVNTLSDSLPGTVGIHGDSPQFFLYSAFMTLVIIILLHVFWGI VFFDG CEKKKWGILLIVLLTHLLVSAQTFISSYYGINLASAFIILVLMGTWAFLAAGGSCR SLKL CLLCQDKNFLLYNQRSR
Source	in vitro E.coli expression system
Target Names	APH1B
Protein Names	Recommended name: Gamma-secretase subunit APH-1B Short name= APH-1b Alternative name(s): Aph-1beta Presenilin-stabilization factor-like
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.