



Recombinant Human Protein FAM57A (FAM57A)

Product Code	CSB-CF844714HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8TBR7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLLTLAGGALFFPGLFALCTWALRRSQPGWSRTDCVMISTRVSSVHAVLATG SGIVIIR SCDDVITGRHWLAREYVWFLIPYMIYDSYAMYLCEWCRTDQNRAPSLTLRN FLSRNRLM ITHHAVILFVLVPVAQRLRGDLGDFVGCIFTAELSTPFVSLGRVLIQLKQQHTLL YKVN GILTATFLSCRILLFPFMYWSYGRQQGLSLLQVPFSIPFYCNVANAFLVAPQIY WFCLL CRKAVRLFDTPQAKKDG
Source	in vitro E.coli expression system
Target Names	FAM57A
Protein Names	Recommended name: Protein FAM57A Alternative name(s): Protein CT120
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.