



Recombinant Escherichia coli UPF0014 inner membrane protein ybbM (ybbM)

Product Code	CSB-CF305194ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77307
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MNSHNITNESLALALMLVVVAILISHKEKLALEKDILWSVGRRAIQLIIVGYVLKYIF SV DDASLTLLMVLFCFNAAWNAQKRSKYIAKAFISSFIAITVGAGITLAVLILSGSIEF IP MQVIPIAGMIAGNAMVAVGLCYNNLGQRVISEQQQIQEKLSLGATPKQASAILIR DSIRA ALIPTVDSAKTVGLVSLPGMMSGLIFAGIDPVKAIKYQIMVTFMLLSTASLSTIIAC YLT YRKFYNSRHQLVVTQLKKK
Source	in vitro E.coli expression system
Target Names	fetB
Protein Names	Recommended name: UPF0014 inner membrane protein ybbM
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.