



Recombinant Sheep Type-2 angiotensin II receptor (AGTR2)

Product Code	CSB-CF631736SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q28929
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	PVLYYIIWGVGFLVNTIVVTLFCCQKGPCKVSSSIYIFNLAVADLLLLATLPLWATY YSHR YDWIFGPVMCKVFGSFLTLNMFASIFFITCMSVDRIYQSVIYPFLSQRNPWQA SYIVPLG WCMACLSLPTFYFRDVRTIEYLGVNACIMAFPPEKYAQWSAGIALMKNILGFII PLIFI ATCYFGIRKHLKTNISYGNRITRDQVLKMAAAVVLAFIICWLPFHVLTFLDALA WMGVI NSCEVIAVIDLALPFAILLG
Source	in vitro E.coli expression system
Target Names	AGTR2
Protein Names	Recommended name: Type-2 angiotensin II receptor Alternative name(s): Angiotensin II type-2 receptor Short name= AT2
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.