



# Recombinant Pongo abelii Cytochrome c oxidase subunit 3 (MT-CO3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-CF015074PYX   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P92696  |
| <b>Product Type</b>      | Transmembrane Protein   |
| <b>Immunogen Species</b> | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)   |
| <b>Sequence</b>          | MAHQSHAYHMKPSPWPLTGALSALLTTSGLTMWFHFHSTLLLLTGLLTNALT<br>MYQWWRD<br>VVRESTYQGHHTLPVQKGLRYGMILFITSEVFFFAGFFWAFYHSSLAPTPQLG<br>GHWPTG<br>IIPLNPLEVPLLNTSVLLASGVSITWAHHSLMENNRTQMIQALLITILLGIYFTLLQ<br>ASE<br>YIEAPFTISDGIYGSTFFMATGFHGLHVIIGSTFLTVCLARQLLFHFTSKHHFGFE<br>AAAW YWHFVDVVWFLYVSIYWWGS |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>Target Names</b>      | MT-CO3  |
| <b>Protein Names</b>     | Recommended name: Cytochrome c oxidase subunit 3 Alternative name(s):<br>Cytochrome c oxidase polypeptide III   |
| <b>Expression Region</b> | 1-261   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.              |