



Recombinant Actinobacillus pleuropneumoniae serotype 7 Prolipoprotein diacylglyceryl transferase (Igt)

Product Code	CSB-CF463109AXI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3H2W4
Product Type	Transmembrane Protein
Immunogen Species	Actinobacillus pleuropneumoniae serotype 7 (strain AP76)
Sequence	MNEQFIQFPQIDPIIFSIGPIALRWYGLMYLIGFGFAYWLGMRRAKNSNGVWTT EQVDQL IYTCFWGVILGGRIGDVFFYNFDRLLQDPMFLFRIWEGGMSFHGGLIGVIVAMI WVSFRQ KRSFWNTADFIAPLIPFGLGMGRIGNFINDELWGRITDVPWAVLFPSSGGYLPRH PSQLYE FFLEGVVLFFILNWFIIKKPRPAGSVAGLFLIGYGVFRFLVEYVRDIDPNVNTVDD LITRG QLLSLPMIIGGLAIMIWAYS RKKA
Source	in vitro E.coli expression system
Target Names	Igt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.