



Recombinant Rhizobium loti Undecaprenyl-diphosphatase 2 (uppP2)

| | |
|--------------------------|---|
| Product Code | CSB-CF859878RCU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q98NJ1 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Rhizobium loti (strain MAFF303099) (Mesorhizobium loti) |
| Sequence | MQGITELLPISSTAHMRIVPALLGWQDPGSAFSAAMQLAALAAVISYFWGDVR DLLFGSL DALTRRDFSDRHFRLASWIVLATIPIVIAGVALSGVLNACNSPLRSLTVIGWSCIA MAIL LALAEIFARHKRTIAEASLADALLVGVAQIGALIPGVSRSGSTLTAALGLGFKRA EAARF SFLGLPAIALAGLKELWELHKVHLDAGWSVLATGLVVASISAFFAIWGLMRV LERFSA WPFVIYRGLLGVVLLLGLAMGWLA |
| Source | in vitro E.coli expression system |
| Target Names | uppP2 |
| Protein Names | Recommended name: Undecaprenyl-diphosphatase 2 EC= 3.6.1.27 Alternative name(s): Bacitracin resistance protein 2 Undecaprenyl pyrophosphate phosphatase 2 |
| Expression Region | 1-264 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |