



Recombinant Pasteurella multocida Prolipoprotein diacylglyceryl transferase (Igt)

Product Code	CSB-CF880385ESG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9CPG2
Product Type	Transmembrane Protein
Immunogen Species	Pasteurella multocida (strain Pm70)
Sequence	<pre>MQSNYFSFPQFDPVIFEIGPIGLRWYGLMYLLGFLFARWLAVKLANQTGSGWT TDQVDSL LFNGFMGVFLGGRIGYVLFYQFDYFLQDPAYLFRVWEGGMSFHGGLIGVILAM LITAKIQ KRGFWQTADFVAPLIPFGLGMGRIGNFINDELWGRVTDVPWAVLFPSSGGYLP RHPSQLYE AVLEGIVLFFILNWIYIKKPRPIGATAGLFLGYGIFRFIVEFFREPDAQLGLYFGQ HISM GQILSTPMILIGAVIMLVAYQSAGKKREIL</pre>
Source	in vitro E.coli expression system
Target Names	Igt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.