



Recombinant *Campylobacter jejuni* subsp. *jejuni* serotype O:6 Prolipoprotein diacylglyceryl transferase (Igt)

Product Code	CSB-CF425520CWG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8FKJ4
Product Type	Transmembrane Protein
Immunogen Species	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> serotype O:6 (strain 81116 / NCTC 11828)
Sequence	MEFWQHIYSNFNVIAFSIFGLKVHWYGYIMYVIALLLALLLAKFFVRKFQLDINEKH LDSY FIWVEIGVILGARLGYILIYDANTMYIYTHPWQIFNPYINGEFVVGIRGMSYHGAIG FLI ATLLFCKKYKTNPWIFLDLVALSVPLAYVFGFRIGNFLNQELFGRITNVPWGIYVD GVLRH PSQLYEAFLEGIVVFIIVYLARFKQSFQGELILVYAGAYSLARFICEFYREPDFGI GFVL WGMSMGQILSFIMFITALLVYICIKFKKVN
Source	in vitro E.coli expression system
Target Names	Igt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.