



# Recombinant Escherichia coli O6:K15:H31 Ferrous iron permease EfeU (efeU)

<b>Product Code</b>	CSB-CF610694EGY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q0TJ50
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Sequence</b>	MFVPPFLIMLREGLEAALIVSLIASYLKRTQRGRWIGVMWIGVLLAAALCLGLGIFI NETT GEFPQKEQELFEGIVAVIAVVILTWMVFWMRKVS RNVKVQLEQAVDSALQRG NHHGWALV MMVFFAVAREGLESVFFLLAAFQQDVGIWPPPLGAMLGLATAVVLGFLLYWGGI RLNLGAF FKWTS LFILFVAAGLAAGAIRAFHEAGLWNHFQEIAFDMSAVLSTHSLFGTLME GIFGYQ EAPSVSEVAVWFIYLIPALVAFVLP PRAGATASRSV
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	efeU
<b>Protein Names</b>	Recommended name: Ferrous iron permease EfeU Alternative name(s): Fe(2+) ion permease EfeU Ferrous iron uptake protein
<b>Expression Region</b>	1-276
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.