



Recombinant Human Putative olfactory receptor 10J6 (OR10J6P)

Product Code	CSB-CF818809HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGY7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRRKNLTEVTEFVFLGFSRFHKHHITLFVVFLILYTLTVAGNAIMTIICIDRHLHT PMY FFLSMLASSKTVYTLFIIPQMLSSFVTQTQPISLAGCTTQTTFFVTLAINNCFLT VMGY DHYMAICNPLRYRVITSKKVCVQLVCGAFSIGLAMAAVQVTSIFTLPFCHTVVG HFFCDI LPVMKLSCINTTINEIINFVVRLFVILVPMGLVFISYVLIISTVLKIASAEGWKKTF A TC AFHLTVVIVHYGCASIAYLMPKSENSIEQDLLLSVT
Source	in vitro E.coli expression system
Target Names	OR10J6P
Protein Names	Recommended name: Putative olfactory receptor 10J6
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.