



Recombinant NADH-ubiquinone oxidoreductase chain 2 (nd2)

Product Code	CSB-CF015077CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24889
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MIVFISLFTLFLTLLSILTNNVIVVWSIFLLMTVVFILLNKSSKSYTSIFNYFVIQESL G LLFLLCSGGLLQFFIILLKIGVAPLHFVIFNVTNNIFNYGLMWFLTFQKLPFLTILL QIF WLSSVYILLFGLLICVQIFVMKSYKNLLIISSTESFNWIVLGVFFSMFNTFYLFYI YFV LMVLLISKFSKTSGYNFINWETTLVFLNIPFSVSFFVKIFSLSEIFKYDSFFTLFLL FTM FLSVLAFSFWLINLSMKNNEETSNNNKMNFIIFPLMVISII
Source	in vitro E.coli expression system
Target Names	nduo-2
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 2 EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 2
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.