



Recombinant Porphyromonas gingivalis Prolipoprotein diacylglyceryl transferase (lgt)

Product Code	CSB-CF752283EYA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7MW41
Product Type	Transmembrane Protein
Immunogen Species	Porphyromonas gingivalis (strain ATCC BAA-308 / W83)
Sequence	MTLPAFITWDFDPVLFITLFGHPVWYGLLFALGLIILGPWIEKMMWEHEKLDISK WFESLA VYVFGTIVGARLGHVLFYDPAYLANPAKIFVTWEGGLASHGGTIGIIIACWLY SRRVT RKSILWVLDRLAVPTGIVAAMIRLGNLTNSEIFGRPTTLPWGFRFIRSEEYRHLV PNMDM GCHPTQIYEALCYLAVFALCMWLYWKRDAARRYSGLIVGVFLTGIFLSRFIIRI KIVQE PWELKLIESVGLNMGQLLSIPFVLAGIWLIIIRAVKNPITQKLS
Source	in vitro E.coli expression system
Target Names	lgt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.