



Recombinant Burkholderia multivorans 4-hydroxybenzoate octaprenyltransferase (ubiA)

Product Code	CSB-CF433449BPS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9AGE2
Product Type	Transmembrane Protein
Immunogen Species	Burkholderia multivorans (strain ATCC 17616 / 249)
Sequence	MLARFPLYLRLVRMDKPIGSLLLLWPTLNALWIASDGHPRWPLIVIFALGTLLMR SAGCA MNDYADRFDRHVKRTADRPLTSGKIRAWEAVAIAVGLSFVAFLLILPLNTLTK ELSVVA LFVAGTYPFMKRFFAIPQAYLGIAFGFGIPMAFAAVQNTVPPLAWVMLIANVFW SIAYDT EYAMVDRDDDIKIGIRTSALTFGRFDVAAVMACYAVTLGIYVWIGIALGFGAAY WVGWAA AAGCAVYHYTLIKGRERMPCFAAFRHNNWLGGVLFAGIAAHYLMAGS
Source	in vitro E.coli expression system
Target Names	ubiA
Protein Names	Recommended name: 4-hydroxybenzoate octaprenyltransferase EC= 2.5.1.- Alternative name(s): 4-HB polyprenyltransferase
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.