



Recombinant Human Mas-related G-protein coupled receptor member G (MRGPRG)

Product Code	CSB-CF774792HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q86SM5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MFGLFGLWRTFDSVVFYLTIVGLGGPVGNGLVLWNLGFRIKKGPFSIYLLHLA AADFLF LSCRVGFSVAQAALGAQDTLYFVLTFWFAVGLWLLAAFSVERCLSDLFPACY QGCRPRH ASAVLCALVWTPPLPAVPLPANACGLLRNSACPLVCPRYHVASVTWFLVLARV AWTAGVV LFVWVTCCSTRPRPRLYGIVLGALLLFFCGLPSVFYWSLQPLLNFLLPVFSPL ATLLAC VNSSSKPLIYSGLGRQPGKREPLRSVLRRALGEGAELGARGQSLPMGLL
Source	in vitro E.coli expression system
Target Names	MRGPRG
Protein Names	Recommended name: Mas-related G-protein coupled receptor member G Alternative name(s): G-protein coupled receptor 169
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.