



Recombinant Schizosaccharomyces pombe Uncharacterized TLC domain-containing protein C17A2.02c (SPAC17A2.02c)

Product Code	CSB-CF519299SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13752
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MSSAVTPSDQGLLPRGPALAEFVKPFCDRFGTLKLRHMHVILLSALFYQIINIL SPVI SRHLSTHYAKLSKCTRLNWDAAHVSSVQSIVLICLGYTCLKEVNAFPDKLFGYS VVAGDI YALTAGYFVWDLYITVRYVHITGIGFVIHAIAALFVITFSYRPLYMYYGPTYLSWE LSTP FLNIHYFLDKTNRTGSKFQMINGFILIVTFICVRIAWGWFSAYSTAIEILNHINVAP WAL SLFYLAANMSLNCLNLFWVSKMIDAIRRRRAHGEKKSTPLQVTSEYAKKNI
Source	in vitro E.coli expression system
Target Names	SPAC17A2.02c
Protein Names	Recommended name: Uncharacterized TLC domain-containing protein C17A2.02c
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.