



Recombinant Salmonella paratyphi A Prolipoprotein diacylglyceryl transferase (lgt)

Product Code	CSB-CF480417SWT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5BFG7
Product Type	Transmembrane Protein
Immunogen Species	Salmonella paratyphi A (strain AKU_12601)
Sequence	MTSSYLHFPDFDPVIFSIGPVALHWYGLMYLVGFVFAMWLAVRRANRPGSGW TKNEVENL LYAGFLGVFLGGRIGYVLFYNFPLFLDNPLYLFRVWDGGMSFHGGLIGVILVMII FARRT KRSFFQVSDFIAPLIPFGLGAGRLGNFINGELWGRVDPNFRFAMLFPGSRAEDI ALLPSH PQWQPIFDTYGVLPRHPSQLYELALEGVVLFILNLFIRKPRPMGAVSGLFLIGY GAFRI IVEFFRQPDAQFTGAWVQYISMGQILSIPMIIAGAIMMVWAYRRRPQQHVS
Source	in vitro E.coli expression system
Target Names	lgt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.