



Recombinant Sheep Adrenocorticotrophic hormone receptor (MC2R)

Product Code	CSB-CF890597SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9TU77
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	MRHILNLYENINSTARNNSDCPAVILPEEIFFTVSIVGVLENLMVLLAVAKNKSLQ SPMY FFICSLAISDMLGSLYKILENVLMIFRNMGYLEPRGSFESTADDVVDLSFILSLLG SICS LSVIAADRYITIFHALQYHSIVTMHRALVVLTVLWAGCTGSGITIVTFSHHVPTVI AFTA LFPLMLAFILCLYVHMFLARSHARRTSSLPKANMRGAILTLVLLGVFIFCWAPF VLHVL LMTFCPADPYCACYMSLFQVNGVLIMCNAVIDPFIYAFRSPELRVAFKKMVICN W
Source	in vitro E.coli expression system
Target Names	MC2R
Protein Names	Recommended name: Adrenocorticotrophic hormone receptor Short name= ACTH receptor Short name= ACTH-R Alternative name(s): Adrenocorticotropin receptor Melanocortin receptor 2 Short name= MC2-R
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.