



Recombinant *Oryza sativa* subsp. *japonica* Bidirectional sugar transporter SWEET13 (SWEET13)

Product Code	CSB-CF647713OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2QR07
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MAGLSLQHPWAFAGLLGNLISFTTYLAPIPTFYRIYKSKSTEGFQSVYPYVVALF SAMLW IFYALIKSNEALLITINAAGCVIETIYIVMYLAYAPKKAKVFTTKILLLLNVGVFGVIL L LTLLLSHGEQRVVSLGWVCVAFSVSVFVAPLSIIKRVIQSRVSVEYMPFSLSLTLT LSAVV WFLYGLLIKDKYVALPNILGFTFGVVQMGLYVFYMNATPVAGEGKEGKGLAA AEELPVV VNVGKLAAATPDRSTGAVHVHPVPRSCAAEAAAEEPEVLVDIPPPPPRAVEV AAV
Source	in vitro E.coli expression system
Target Names	SWEET13
Protein Names	Recommended name: Bidirectional sugar transporter SWEET13 Short name= OsSWEET13
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.