



Recombinant Mouse Adrenocorticotrophic hormone receptor (Mc2r)

Product Code	CSB-CF723737MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q64326
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MKHIINSYEHTNDTARNNSDCPDVVLPEEIFFTISVIGILENLIVLLAVIKNKNLQS PMY FFICSLAISDMLGSLYKILENILIMFRNMGYLKPRGSFESTADDIIDCMFILSLLGSI FS LSVIAADRYITIFHALQYHSIVTMRRTIITLTIWMFCTGSGITMVIFSHHIPTVLTFT S LFPLMLVFILCLYIHMFLARSHARKISTLPRNTMKGAMTLTILLGVFIFCWAPFV LHVL LMTFCPNNPYCVCYMSLFQVNGMLIMCNAVIDPFIYAFRSPELRDAFKRMLFC NRY
Source	in vitro E.coli expression system
Target Names	Mc2r
Protein Names	Recommended name: Adrenocorticotrophic hormone receptor Short name= ACTH receptor Short name= ACTH-R Alternative name(s): Adrenocorticotropin receptor Melanocortin receptor 2 Short name= MC2-R
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.