



Recombinant Arabidopsis thaliana GDT1-like protein 1, chloroplastic (At1g64150)

Product Code	CSB-CF822006DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q94AX5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	CFQRNDAACYLEKAESEEHDRNLDVLVESSIAHSRREIQRVLMFLAVSGSVAL LGTDPAF AASSIPNVTQSLVTSFGDLGDISSGFASAFLLIFFSELGDKTFFIAALLAARNNSAA TVFV GTFGALGIMTIISVVLGRTFHYVDEVLPFRFGGTDLPIDDIAAVCLLVYFGVSTLL DAVS DEGKADEEQKEAELAVSELSGNGAGIVAAANTIISTFALVFVAEWGDKSFFSTI ALAAAS SPLGVIAGALAGHGAATLLAVLGGSLGNFLSEKAIAYVGGVFLVFAAVTVAEI VT
Source	in vitro E.coli expression system
Target Names	PAM71
Protein Names	Recommended name: GDT1-like protein 1, chloroplastic
Expression Region	74-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.