



Recombinant Human Putative olfactory receptor 5AK3 (OR5AK3P)

Product Code	CSB-CF836774HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH89
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGRGNSTEVTEFHLLGFGVQHEFQHVLFIIVLLLIYVTSIGNIGMILLIKTDSRLQ TPMY FFPQHIAFVDICYTSAITPKMLQSFTENNLIIFRGCVIQFLVYATFATSDCYLLA IMAM DCYVAICKPLRYPMIMSQTVYIQLVAGSYIIGSINASVHTGFTFSLSFCKSNKINH FFCD GLPILALSCSNIDINIILDVVFVGFDLMFTELVIIFSIIYIMVTILKMSSTAGRKKSF S T CASHLTAVTIFYGTLSYMYLQPQSNNSQENMKVASIFYGTVIPMLNPLIYSLRN KEGK
Source	in vitro E.coli expression system
Target Names	OR5AK3P
Protein Names	Recommended name: Putative olfactory receptor 5AK3
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.