



Recombinant Pan troglodytes Taste receptor type 2 member 50 (TAS2R50)

Product Code	CSB-CF731072EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q646C3
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MITFLYIFFSILIMVLFVVLGNFANGFIALVNFIDWVKRKKISSADQILTALAVSRIGL LW TLLLNWYLTVLNPAFYVELRITSYNAWVVTNHFSMWLAASLSIFYLLKIANFSN LIFLH LKRRVRSVILVILLGTLIFLVCHLLVANMDESMWAEYEGNITGKMCLRNTVHLS YLTVT TLWSFIPFTLSLISFLMLICSLCKHLKMKQLHGEGSQDLSTKVHIKALQTLISFLLL CAI FFLFLIISVWSPRRLRNDPVVMVSKAVGNIYLAFDSFILIWRTKCLKHTFLLILCQI RC
Source	in vitro E.coli expression system
Target Names	TAS2R50
Protein Names	Recommended name: Taste receptor type 2 member 50 Short name= T2R50
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.