



Recombinant Mesocricetus auratus Free fatty acid receptor 1 (FFAR1)

Product Code	CSB-CF008604MRG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q76JU8
Product Type	Transmembrane Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Sequence	MALSPQLFFALYVSAFALGFPLNLLAIRGAVARARLRLTPNLVYTLHLACSDLLL AITLP VKAVEALASGAWPLPLPLCPVFLVHFAPLYAGGGFLAALSAGRYLGAAPFPG YQAVRRP RYSWGVCVAIWALVLCHMGLVLGLEAPGGWLNTTSSSLGINTPVNGSPVCLE AWDPNSAR PARLSFSILLFFVPLVITAFICYVGCLRALAHSGLSHKRKLRAAWAAGGAFLTLL CLGPY NASNVASFVNPDLGGSWRKLGLITGSWSVVLNPLVTGYLGASPGRGTVCTTR TQGGTIQK
Source	in vitro E.coli expression system
Target Names	FFAR1
Protein Names	Recommended name: Free fatty acid receptor 1 Alternative name(s): G-protein coupled receptor 40
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.