



Recombinant Human Post-GPI attachment to proteins factor 3 (PGAP3)

Product Code	CSB-CF856935HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96FM1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SQGDREPVYRDCVLQCEEQNCSGGALNHFRSRQPIYMSLAGWTCRDDCKYE CMWVTVGLY LQEGHKVPQFHGKWPF SRFLFFQEPASAVASFLNGLASLVMLCRYRTFVPAS SPMYHTCV AFAWVSLNAFWSTVFHTRDTDLTEKMDYFCASTVILHSIYLCCVRTVGLQHP AVVSAFR ALLLLMLTVHVSYSLSLIRFDYGYNLVANVAIGLVNVVWWLAWCLWNQRRLPHV RKC VVVV LLLQGLSLLELLDFPPLFWVLDAAHAIWHISTIPVHVLFFSFLEDDSLYLLKESEDK FKLD
Source	in vitro E.coli expression system
Target Names	PGAP3
Protein Names	Recommended name: Post-GPI attachment to proteins factor 3 Alternative name(s): COS16 homolog Short name= hCOS16 Gene coamplified with ERBB2 protein PER1-like domain-containing protein 1
Expression Region	21-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.