



Recombinant Pan troglodytes Taste receptor type 2 member 13 (TAS2R13)

Product Code	CSB-CF723757EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q646B6
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MESALPSIFTLVIAEFIIGNLSNGFIVLINCIDWVSKRELSSVDKLLIILAI SRIGLIW EILVSWFLALHSLAIFVSGTGLRIMIFSWIVSNHFNLWLATILSIFYLLKIASFSSPA FL YLKRRVNVKILMILLGTLVFLFLNLIQINMLIKDWLD RYERNTTWNFMSDFETF SVSVR FTMTMFSLTPFTVAFISFLLL VFSLQKHLQKMLNYKGHRDPRTKVHTNALKIVI SFLLL YASFFLSILISWISELYQNTVIYMLCETIGAFYPSSHSFLLILGNAKLRQAFLLVAA KVV AKR
Source	in vitro E.coli expression system
Target Names	TAS2R13
Protein Names	Recommended name: Taste receptor type 2 member 13 Short name= T2R13
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.