



Recombinant Gorilla gorilla gorilla Taste receptor type 2 member 13 (TAS2R13)

Product Code	CSB-CF731063GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q645Z5
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	MKSALPSIFTLVIIAEFIIGNLSNGFIVLINCIDWVSKRELSSVDKLLIILAI SRIGLIW EILVSWFLALHYLAIFVSGTGLRIMIFSWIVSNHFNLWLATILSIFYLLKIASFSSPA FL YLKWRV NKVILLILLGTLVFLFLNLIQINMHIKD WLD RYERNTTWNFMSDFETF SVSVK FTMTMFSLTPFTVAFISFLLLIFSLQKHLQKMLNYKGHRDPKTKVHTNALKIVIS FLLF YASFFLCVLXSWISELYQNTVIYMLCETIGVFYFPSSH S FLLILGNAKLRQAFLLVA AKVW AKR
Source	in vitro E.coli expression system
Target Names	TAS2R13
Protein Names	Recommended name: Taste receptor type 2 member 13 Short name= T2R13
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.