



Recombinant Human Olfactory receptor 4X2 (OR4X2)

Product Code	CSB-CF836751HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGF9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTEFIFLVLSNPQEVQRVCFVIFLFLYTAIVLGNFLIVLTVMTSRSLGSPMYFFLS YLSF MEICYSSATAPKLISDLLAERKVISWWGCMAQLFFLHFFGGTEIFLLTVMAYDH YVAICK PLSYTTIMNWQVCTVLVGIWVGGFMHSFAQILLIFHLLFCGPNVINHYFCDLV PLLKLA CSDTFLIGLLIVANGGTLVISFGVLLASYMVILLHLRTWSSEGWCKALSTCGSH FAVVI LFFGPCVFNSLRPSTTLPIDKMVAVFYTVITAILNPVIYSLRNAEMRKAMKRLWI RTLRL NEK
Source	in vitro E.coli expression system
Target Names	OR4X2
Protein Names	Recommended name: Olfactory receptor 4X2 Alternative name(s): Olfactory receptor OR11-105
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.