



Recombinant Human Olfactory receptor 10A2 (OR10A2)

Product Code	CSB-CF884445HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H208
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSFSSLPTEIQSLLFLFTLIYLVTLMGNCLLILVTLADPMLHSPMYFFLRNLSFLE IGF NLVIVPKMLGTLAQQDTTISFLGCATQMYFFFFFFGVAECFLATMAYDRYVAICS PLHYP VIMNQRTRAKLAAASWFPGFVATVQTTWLFSFPFCGTNKVNHHFFCDSPPVL RLVCADTA LFEIYAIVGTILVVMIPCLLILCSYTHIAAAILKIPSAKGKNAFSTCSSHLLVVSFLF YI SLSLTYFRPKSNNSPEGKLLSLSYTVMTPLNPIIYSLRNNEVKNALSRTVSK ALALRN CIP
Source	in vitro E.coli expression system
Target Names	OR10A2
Protein Names	Recommended name: Olfactory receptor 10A2 Alternative name(s): HP4 Olfactory receptor OR11-86
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.