



Recombinant Human Putative olfactory receptor 52A4 (OR52A4)

Product Code	CSB-CF016818HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NMU1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MALPITNGTLFMPFVLTFIGIPGFESVQCWIGIPFCATYVIALIGNSLLLIIKSEPSL H EPMYIFLATLGATDISLSTSIVPKMLDIFWFHLPEIYFDAFLFQMWLIHTFQGIES GVLL AMALDRCVAICYPLRRAIVFTRQLVTYIVVGVTLRPAILVIPCLLLIKCHLKLYRTK LIY HTYCERVALVKLATEDVYINKVYGILGAFIVGGLDFIFITLSYIQIFITVFHLPLKEA RL KVFNTCIPHIYVFFQFYLLAFFIFYSQIWILYPIICTYHLVQSLPTGPTIPQPLYLW VK DQTH
Source	in vitro E.coli expression system
Target Names	OR52A4P
Protein Names	Recommended name: Putative olfactory receptor 52A4
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.