



Recombinant Mouse Probable G-protein coupled receptor 141 (Gpr141)

Product Code	CSB-CF749267MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TQP0
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MDGYNTSENSSCDPILAHHLTSIYFIVLIGGLVGLISILFLLVKMNSRSVTTMAVIN LVV VHGVFLLTVPFRLAYLIKGTWTFGLPFCKFVSAMLHIHMYLTFLFYVVILVIRYLIF FKR RDKVEFYRKLHAVAASSAMWLLVIVIVVPLVVSQYGNSEEYNEQQCFRFHKEL GHDSVRV INYMIVIVVIAVALILLGFQVFITLSMVRKFRHSLLSHQEFWAQLKNLFFIGIIIICFL P YQFFRIYYLYVVAHSKSKKNKVAFYNEILLSTTAISCCDLLLFVFGGSHWVKQKI VDMWN CLLCH
Source	in vitro E.coli expression system
Target Names	Gpr141
Protein Names	Recommended name: Probable G-protein coupled receptor 141 Alternative name(s): G-protein coupled receptor PGR13
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.