



Recombinant Human Olfactory receptor 9G9 (OR9G9)

Product Code	CSB-CF017213HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C7N8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MQRSNHTVTEFILLGFTTDPGMQLGLFVVFLGVYSLTVVGNSTLIVLICNDSHL HTPMYF VVGNSLFLDLWYSSVYTPKILVICISEDKISFAGCLCQFFFSAGLAYSECCLA AMAYD RYVAISKPLLYAQAMSIKLCALLVAVSYCGGFINSSIITKKTFSFNFCENIIDDF CDL LPLVKLACGEKGCYKFLMYFLLASNVICPAVLILASYLFIITSVLRISSSQGRLKAF STC SSHLTSVTLYYGSILYIYALPRSSYSFDMDKIVSTFYTEVLPMLNPMIYSLRNKD VKEAL KKLLP
Source	in vitro E.coli expression system
Target Names	OR9G9
Protein Names	Recommended name: Olfactory receptor 9G9
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.