



Recombinant Papio hamadryas Taste receptor type 2 member 10 (TAS2R10)

Product Code	CSB-CF723761ERJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q646F5
Product Type	Transmembrane Protein
Immunogen Species	Papio hamadryas (Hamadryas baboon)
Sequence	MLSVVEGILILVVISVFGVVGNGFIGLVNCIDCAKNKLSTIGFILTGLAISRIFLI WI IITDGFIIQIFSPDVYASGNLIEYISYFWVITNQSSIWVATSLSIFYFLKIANFSNYIFL W LKSRINRVLPLLMGFLNISCLLNFAIVKILNDLKMKNNTVWRLNMYKSEYFIKQL LLNL GVIFFTLSLITSVLLIISLWRHNRQMMSNVTGLRDSITEAHVKAMKVLISFIILFIL YF IGIAIEISYFTVPENKLLLIFGMTTITAIYPWGHFSFILILGNSKLKQASLRVLQQLKC CEE RKNLRAT
Source	in vitro E.coli expression system
Target Names	TAS2R10
Protein Names	Recommended name: Taste receptor type 2 member 10 Short name= T2R10
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.