



Recombinant Human Taste receptor type 2 member 41 (TAS2R41)

Product Code	CSB-CF023159HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59536
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MQAALTAFFVLLFSLLSLLGIAANGFIVLVLGREWLRYGRLLPLDMILISLGASRF CLQL VGTVHNFYSSAQKVEYSGGLGRQFFHLHWHFLNSATFWFCSWLSVLFCVKIA NITHSTFL WLKWRFPGWVPWLLLGSVLISFIITLLFFWVNYVPVYQEFLIRKFSGNMTYKWN TRIETYY FPSLKLVIWSIPFSVFLVSIMLLINSLRRHTQRMQHNHSLQDPSTQAHTRALK SLISFL ILYALSFLSLIIDAAKFISMQNDFYWPWQIAVYLCISVHPFILIFSNLKLRVFSQLL LL ARGFWVA
Source	in vitro E.coli expression system
Target Names	TAS2R41
Protein Names	Recommended name: Taste receptor type 2 member 41 Short name= T2R41 Alternative name(s): Taste receptor type 2 member 59 Short name= T2R59
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.