



Recombinant Human Olfactory receptor 4D2 (OR4D2)

Product Code	CSB-CF016700HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58180
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	METGNLTWVSDFVFLGLSQTRELQRFLFLMFLFVYITTVMGNILIIITVTSDSQLH TPMY FLLRNLAVIDLCFSSVTAPKMLVDLLSEKKTISYQGCMGQIFFFHFLGGAMVFF LSVMAF DRLIAISRPLRYVTVMNTQLWVGLVVATWVGGFVHSIVQLALMLPLPFCGPNIL DNFYCD VPQVLRRLACTDTSLLLEFLKISNSGLLDVWVFFLLMSYLFILVMLRSHPGEARRK AASTC TTHIIVVSMIFVPSIYLYARPFTPFPMDKLVSIGHTVMTPLNPMIYTLRNQDMQ AAVRR LGRHRLV
Source	in vitro E.coli expression system
Target Names	OR4D2
Protein Names	Recommended name: Olfactory receptor 4D2 Alternative name(s): B-lymphocyte membrane protein BC2009 Olfactory receptor OR17-24
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.