



Recombinant Gamma-secretase subunit aph-1 (aph-1)

Product Code	CSB-CF523691CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O45876
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MGYLLTIACYIASFSPSIALFCSFIAHDPVRIILFFLGSSFFWLVSLLFSSLAWLGLS TVL PDTFLLSLTVCIIAQELSRVAYFMLLKKAQRGLNKITRQGQISVAPGVSDLHNAR HMLAL VCGLGMGVISALFYTMNAFAIFSGPGTIGLPNALKTGEIDTNRAGKYLPLCYTLS AILLT LFHVTWTIMVWDSCHKIGRIPSAFVPGAAAVVSHLLVTFLSSLNSRGFHVLFVA VQFLIL LICIAYCNVIMGGTISSFVNGIGQSITDAVTLKQVRTLIEERKLRTQRQSVDPDEP MTERA GTSNTVNA
Source	in vitro E.coli expression system
Target Names	aph-1
Protein Names	Recommended name: Gamma-secretase subunit aph-1 Alternative name(s): Anterior-pharynx-defective protein 1 Presenilin enhancer protein 1
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.