



Recombinant Macaca mulatta Taste receptor type 2 member 46 (TAS2R46)

Product Code	CSB-CF734117MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q645T4
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	MITFLSITFSILVGVIFVIGNFANGFIALVNSIEWVKRQKISFADQILTGLAVSRVGL LW VLLLHLYATEFNLAFY SVEVRITAYNVWIVTNHFSNWLSTLSM FYLLKIATFSN LIFLH LKRKVKSVILVTLLGPLLFLVCHL FVMNMNHIVWRKEYEGNITWRIKLRSAMYL SNVTVT MLANLIPLTLT LMSFLLLICSLCKHLKMKQVHGKGSQDPSTKVHIKALQTVTSFL LLCAI YFLSMILSVWNFELEK KPVFMFCQAVIFSYPSTHPLILIWGNK KKLKQIFLSVLWN VRYWV KGQKPSSP
Source	in vitro E.coli expression system
Target Names	TAS2R46
Protein Names	Recommended name: Taste receptor type 2 member 46 Short name= T2R46
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.