



# Recombinant Human Olfactory receptor 2T6 (OR2T6)

|                          |                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF839915HU                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                          |
| <b>Uniprot No.</b>       | Q8NHC8                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product Type</b>      | Transmembrane Protein                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                       |
| <b>Sequence</b>          | MNENNETLTRGFTLMGLFTHNKCSGFFFGVICAVFFMAMIANGVMIFLINIDPHL<br>HTPMY<br>FLLSHLSVIDTLYISTIVPKMLVDYLMGEGTISFIAC TAQCFLYMGFMGAEFFLLG<br>LMAY<br>DRYVAICNPLRYPVLISWRVCWMILASSWFGGALDSFLLTPITMSLPFCASHQI<br>NHFFCE<br>APTMLRLACGDKTTYETVMYVCCVAMLLIPFSVVTASYTRILITVHQMTSAEGR<br>KKAFAT<br>CSSHMMVVTLFYGAALYTYTLPQSYHTPIKDKVFSAFY TILTPLLNPLIYSLRNR<br>DVMGA LKRVARC |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                                                                                          |
| <b>Target Names</b>      | OR2T6                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Names</b>     | Recommended name: Olfactory receptor 2T6 Alternative name(s): OST703<br>Olfactory receptor 2T9                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-308                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                        |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                              |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                         |