



Recombinant Mouse Cysteinyl leukotriene receptor 2 (Cysltr2)

Product Code	CSB-CF006466MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q920A1
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEVTGTPSSYSNRNCTIENFKKEFYPIIYLIIFFWGALGNGFSIYVFLQTCKKSTS VNVF MLNLATSDFLFISTLPFRADYYFRGSNWIFGDLACRVMSYSLYVNMYSIYFLT VLSVVR FLATVHPFRMFHVTSVRSAWILCGIIVWFIMASSALLLVNGQEEKDNIISCLELS PQKFK SLLIMNHIAVAVGFLLPFLTLTVCYLLIIRILLKAEIPESGPRAAHRKALTTIVIAMIT F LLCFLPYHALRTLHLVTWDKDSCGDVLHKATVITLTMAAANSCFNPFLYYFAGE NFKARL RAIFSKVHL
Source	in vitro E.coli expression system
Target Names	Cysltr2
Protein Names	Recommended name: Cysteinyl leukotriene receptor 2 Short name= CysLTR2
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.