



# Recombinant Pan paniscus Taste receptor type 2 member 20 (TAS2R20)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-CF734130EQU   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | Q646E3  |
| <b>Product Type</b>      | Transmembrane Protein   |
| <b>Immunogen Species</b> | Pan paniscus (Pygmy chimpanzee) (Bonobo)  |
| <b>Sequence</b>          | MMSFLHIVFSILVVVAFILGNFANGFIALINFIWVVKRQKISSADQIIAALAVSRVGLLW<br>VILLHWYSTVLNPTSSNLKVIIFISNAWAVTNHFSIWLATLSLIFYLLKIVNFSRLIFHH<br>LKRKAKSVVLVIVLGLSFLVCHLVMKNTYINWVTEECEGNVTWKIKLRNAMHL<br>SNLTVA<br>MLANLIPFTLTLISFLLLIYSLCKHLKMKMLHGKGSQDPSTKIHIALQTVTSFLILLAI<br>YFLCLITSFWNSKMRPKEIVLMLCQAFGIIYPSFHSFILIWGNKTLKQTFLSVLWQVTCW AKGQNGQSTP |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>Target Names</b>      | TAS2R20   |
| <b>Protein Names</b>     | Recommended name: Taste receptor type 2 member 20 Alternative name(s): Taste receptor type 2 member 49 Short name= T2R49  |
| <b>Expression Region</b> | 1-309   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |