



Recombinant Rat Mas-related G-protein coupled receptor member E (Mrgpre)

Product Code	CSB-CF742495RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7TN40
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSLRVHTHSPSTQGDMAFNLTILSLTELLSLGGLLGNGVALWLLNQNQVYRNPF SIYLLDV ACADLIFLCCHMVAIPELLQDQLNFPEFVHISLIMLRFFCYIVGLSLLVAISTEQC LAT LFPSGYLCRRPRYLTTVCFAFIWVLCLLLDLLSGACTQFFGAPSYHLCGMLW LVVAVLL AALCCTMCVTSLLLLLRVERGPERHQPRGFPTLVLLVILLFLFCGLPFGIFWLSK NLSWH TPLYFYHFSFFMASVHSAAKPAIYFFLGSTPGQRFQEPLRLVLQRALGDEAEL GAVREAS QGGLVDMTV
Source	in vitro E.coli expression system
Target Names	Mrgpre
Protein Names	Recommended name: Mas-related G-protein coupled receptor member E
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.