



Recombinant Mouse Protein lifeguard 3 (Tmbim1)

Product Code	CSB-CF805910MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BJZ3
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSNPSAPPPYEDHNPLYPGSPPPGGYGQPSVLPGGYPAYPAYPQPGYGHPA GYPQPVPV HPMPMNYGHDYNEEERAGSDSFRPGEWDDRKVRHSFIQKVYCIISVQLLITVA IIAIFTF VEPVGKYVRNNVAVYYVSYAVFLVITYLTLACCQGPRRRFPWDIILLTIFTLALGF VTGTI SSMYENKAVIIAMIITAVVSISVTIFCFQTKVDFTSCTGLFCVLGIVLMVTGIVTSIV LI FKYIYWLHMVYAALGAICFTLFLAYDTQLVLGNRKHTISPEDYITGALQIYTDIVYI FTF VLQLVGSRD
Source	in vitro E.coli expression system
Target Names	Tmbim1
Protein Names	Recommended name: Protein lifeguard 3 Alternative name(s): Responsive to centrifugal force and shear stress gene 1 protein Short name= Protein RECS1 Transmembrane BAX inhibitor motif-containing protein 1
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.