



Recombinant Human Olfactory receptor 4C12 (OR4C12)

Product Code	CSB-CF847215HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96R67
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEKKKNVTEFILIGLTQNPIMEKVTFVFLVLYMITLSGNLLIVVTITTSQALSSPM YFF LTHLSLIDTVYSSSSAPKLIVDSFQEKKIISFNGCMAQAYAEHIFGATEIILLTVMA CDC YVAICKPLNYTTIMSHSLCILLVAVAWVGGFLHATIQLFTVWLPFCGPNVIGHF MCDLY PLLKLVCIDHTLGLFVAVNSGFICLLNFLILVVSYVILRSLKNNSLEGRCKALST CIS HIIVVLFVPCIFVYLRSVTTLPIDKAVAVFYTMVVPMLNPVYTLRNAEVKSAI RKLW RKKVTSDND
Source	in vitro E.coli expression system
Target Names	OR4C12
Protein Names	Recommended name: Olfactory receptor 4C12 Alternative name(s): Olfactory receptor OR11-259
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.