



Recombinant Leontopithecus chrysomelas Melanocyte-stimulating hormone receptor (MC1R)

Product Code	CSB-CF801135LLN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q864I1
Product Type	Transmembrane Protein
Immunogen Species	Leontopithecus chrysomelas (Golden-headed lion tamarin)
Sequence	MPMQGAQRKLLGSLNSTPTATSNLGLAANRTGAPCLELPIPDGLFLSLGLVSL VENLVV AAIKRNRLHSSMYCFICCLALSDLLVSGSNMLEAGVLATRASVVQQLHNTIDV LTCSSM LCSLCFLGAIAVDRIYISIFYALRYHSIMTLPRARAAIWWASVLSSTLFITYYD HAAV LLCLVVFFLAMLVLMVLYVHMLAWACQHAQGIIRLHKRQPPAHKGFGLRGAA TLTILLG IFFLCWGPFFLRRLTLVVFQHLTCNCIFKNFKVFLTLICNTIIDPLIYAFRSQELR RT LKEVLGRGRW
Source	in vitro E.coli expression system
Target Names	MC1R
Protein Names	Recommended name: Melanocyte-stimulating hormone receptor Short name= MSH-R Alternative name(s): Melanocortin receptor 1 Short name= MC1-R
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.