



Recombinant Mouse Olfactory receptor 472 (Olf472)

Product Code	CSB-CF848702MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VG43
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEAENHTTVAELIILGLTEDPKLCIVFFVIFLGVYIITLVGNISIITLIRISSQLHTPMY LFLSHLAFVDIVFSTSVSVIMLMELLGHGLVLSVATCAAQLCMTVSFGSAECFLL AAMAY DRYVAICSPLLYSTLMSSRVCFLLLGISYVGGFVNGWTFGTGCVLSLSFCGPTQI NHFFCD FSPLLKVSCSDVSIIGIIPSISSSGSIIVTVFVIAVSYIYLITILKMRSTEGRHKAFST CTSHLTAVTLFYGTITVIYVMPKSSSYSTEQNKVISLFYTVVIPMLNPLIYSLRNRD VKDA LRKAIVRVYS
Source	in vitro E.coli expression system
Target Names	Olf472
Protein Names	Recommended name: Olfactory receptor 472 Alternative name(s): Olfactory receptor 204-5
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.