



Recombinant Human Olfactory receptor 8D2 (OR8D2)

Product Code	CSB-CF867107HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9GZM6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MATSNHSSGAEFILAGLTQRPELQLPLFLLFLGIYVVTVVGNLGMIFLIALSSQLY PPVY YFLSHLSFIDLICYSSVITPKMLVNFVPEENIISFLECITQLYFFLIFVIAEGYLLTAM EY DRYVAICRPLLYNIVMSHRVCSIMMAVVYSLGFLWATVHTTRMSVLSFCRSHT VSHYFCD ILPLTLSCSSTHINEILLFIIGGVNTLATTLAVLISYAFIFSSILGIHSTEGQSKAFG T CSSHLLAVGIFFGSITFMFKPPSSTTMEKEKVSSVFYITIIPMLNPLIYSLRNKD VKNA LKKMTRGRQSS
Source	in vitro E.coli expression system
Target Names	OR8D2
Protein Names	Recommended name: Olfactory receptor 8D2 Alternative name(s): Olfactory receptor OR11-303 Olfactory receptor-like protein JCG2
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.